

STYL'One Evolution

R&D, SCALE-UP AND PRODUCTION SUPPORT

+
40%
Output

+
60%
Speed

NEW EVOLUTION OF STYL'ONE TABLETING INSTRUMENT THE ALL-IN-ONE MULTIFUNCTIONAL TOOL FROM EARLY FORMULATION TO SCALE-UP AT PRODUCTION SPEED


- Single layer tablets production upgradable to multi-layer and tab-in-tab
- Target thickness or target density mode
- Rotary Press Mimicking – typically up to 400 000 tab/ hour (Dwell time <2 ms)
- Dry granulation with roll compaction hybrid modeling
- Full library of rotary tablet presses and roll compactors
- Press Coating Device (tab-in-tab)
- Compact Design, door size, easy to move
- Clean & easy tooling (B&D) change
- Extensive data analysis and automatic studies with only a minimal amount of powder
- Interchangeable die table
- Premium instrumentation



STYL'ONE EVOLUTION | TECHNICAL DATA



Extensive vacuum system




Interchangeable die table



Multilayer | Tab-in-Tab

Number of press stations	1
Tool type	B+D (EU/ TSM), oversized, core rods
Die	BBS, BB, B, D, oversized
Maximum tablet diameter (mm)	21 (B) – 25 (D) – 40 (oversized)
Maximum output (tablets/ hour)	1 750
Feed shoes (up to 3)	Hand feed, gravity, paddle forced feeder, core feeding
Maximum die filling (mm)	21 (B) – 23 (D) – 34 (oversized)
Compression mode	Force or displacement driven
Compression zone	Adjustable 1 to 12 mm insertion depth.
Max. pre-compression force (kN)	50 kN standard *
Max. compression force (kN)	50 kN standard *
Punch movement	Selectable 1 punch only or upper and lower
Dwell time (ms)	<2 ms can be extended up to >2 sec
R&D accessories	Instrumented die, heated die, enhanced operator safety, full containment, automated studies and more
Production accessories	GMP Package, weight regulation, external lubrication and more
Power (kW)	2 kW max during operation
Standard Voltage	400/ 480 V (± 10%) – 50/60 Hz
Weight (kg)	1090 kg (2 403 lbs)

* if other maximum is required, please consult us



MEDELPHARM SAS

12, rue des Petites Combes
F – 01700 Beynost (Lyon), France
Phone +33 478 97 62 10 | Fax +33 478 88 18 65
contact@medelpharm.com | www.medelpharm.com